


What's in the box?




Eye Socket




Face Mask



Controller



Charger Cable



Velcro Strap


Important Information

Make sure to review the information below before using the Nanoleaf LED Face Mask.



Warning


This is not a medical device and cannot be used for professional treatment or as a substitute for professional medical advice. The actual effect of use varies from person to person.





1. **Do NOT** use the device while it's charging. If the device becomes excessively hot, emits smoke, or produces an unusual odor during charging, please immediately disconnect the charger and contact customer service.



2. Avoid folding, squeezing, or pressing down on the mask to prevent damage.



3. It's important to use eye covers to protect your eyes during treatment. Avoid looking directly at the LED lights for safety and comfort while wearing the mask.




4. When cleaning the mask, do not soak it in water. Simply wipe it with a damp cloth or using 70% alcohol.

Sensitivity Test

It is important to conduct a sensitivity test to ensure the mask is suitable for your skin.

1. Follow the "How to use Nanoleaf LED Face Mask" steps outlined below.
2. Let the device run for the full 10-minute treatment cycle. It will automatically switch off once the session is complete.
3. If at any point your skin feels hot or uncomfortable, stop using the device immediately.
4. After the treatment, wait at least 6 hours. If you notice any redness or itchiness, the device may not be suitable for your skin. In this case, you can return it within 30 days for a full refund.

Before Using the Mask



Ensure your mask is fully charged before its first use. It takes approximately 3.5 hours to reach a full charge. The mask operates for a total of 80 minutes on a full charge before requiring recharging.


If the charging symbol on the controller's screen is flashing and the controller beeps continuously for 2 seconds, the battery is depleted and it's time to charge.

How to Use the Nanoleaf LED Face Mask

1

Attach the Straps

Insert straps and tighten with velcro.




2

Install the Eye Sockets

Gently insert the eye sockets from inside the mask.

Note

This step is important to protect your eyes during treatment, ensuring safety and comfort while using the mask.



3

Power On

Press and hold the power button (⏻) for 2 seconds until you hear an audible sound.

4

Pause Timer

After turning on the mask, the timer will automatically start counting down from the default setting of 10 minutes. Single press the power button (⏻) once to pause the timer.

5

Choose Mode

To select your preferred treatment mode (color), press the arrow button (➡). The colors cycle in the following order: Red > Green > Blue > Cyan > Yellow > Purple > White.

6

Choose Level of Intensity

Double press the arrow button (➡) to cycle through the intensity levels which range from 1 (lowest) to 5 (highest).

7

Choose Time

The default timer is set to 10 minutes. You can increase the time by pressing the timer button (⌚), with each press adding an additional minute. The timer can be set between a minimum of 1 minute and a maximum of 30 minutes. **To select a time from 1 to 9 minutes, simply cycle past 30 minutes, and the timer will reset to 1 minute.**

The face mask will automatically turn off once the session is complete.

8

Put the mask on & enjoy your relaxation time!



Nanoleaf LED Face Mask Controls



controller



Power Button

By default, your device is set to turn off after 10 minutes. Press the power button once more when on to pause the timer.

Hold for 2s to Turn On / Off



Mode Selection

Each light mode is combined with Near Infrared Light (NIR) to enhance the treatment effects. The screen will display the corresponding light mode.

Press once to switch between Red, Green, Blue, Purple, Cyan, Yellow, & White light therapy.

Press twice quickly to cycle between 5 different levels of intensity. The screen will display the corresponding light intensity.



Timer

Your device will automatically turn off when your timer is up.

Press once to add 1 minute to your timer.

Long press to add more minutes quickly up to a max of 30 minutes.

Light Therapy Modes & Treatment Time

1

Light Therapy Modes

Each light mode is combined with Near Infrared Light (NIR) that has 850nm penetration to enhance the treatment effects.



Mode R: Red Light (640nm + 850nm): Increases cell activity and promotes collagen formation to improve skin elasticity, and reduce the appearance of wrinkles



Mode G: Green Light (520nm + 850nm): Enhances the stimulation of cell energy to help even out skin tone, reduce redness and hyperpigmentation.



Mode B: Blue Light (460nm + 850nm): Antibacterial and anti-inflammatory to help reduce acne breakouts and repair the skin.



Mode C: Cyan Light (460nm + 525nm + 850nm): Soothes and calms stressed out skin. Promotes healing and treats inflammation.



Mode Y: Yellow Light (640nm + 525nm + 850nm): Promotes tissue healing to help treat rosacea, reduce redness and improve overall skin tone.



Mode P: Purple Light (640nm + 460nm + 850nm): Multi-frequency treatment that targets acne and acne scars.



Mode L: White Light (640nm + 525nm + 460nm + 850nm): Penetrates the skin deeper to accelerate skin metabolism and enhance nutrient absorption.

2

Recommended Treatment Time for Each Mode

Optimize your skincare results by following the suggested durations for each mode, carefully calibrated for safety and effectiveness.

Aim for **5-10 minutes per mode**:

- **5 minutes** provides a moderate dosage (30J/cm²), ideal for most skincare routines.
- **10 minutes** ensures deeper penetration at full intensity for enhanced results.

The recommended optimal times for each mode are:

 **Mode R** 6mins

 **Mode G** 7mins

 **Mode B** 6mins

 **Mode C** 5mins

 **Mode Y** 6mins

 **Mode P** 6mins

 **Mode L** 6mins

Choose the duration that fits your skincare goals and enjoy effective, tailored treatments!

FAQs

1

How does LED Therapy work?

LED Therapy—or photobiomodulation— works by delivering different wavelenths of light to the skin to stimulate cells to work more efficiently.

The light energy boosts the production of cellular energy, which helps cells repair damage, produce collagen, and improve circulation. This process supports the skin's natural healing and rejuvenation, addressing issues like fine lines, acne, and dullness in a non-invasive way.

The Nanoleaf LED Light Therapy Face Mask comes with 7 different light therapy modes.

2

What is Near Infrared Light (NIR)?

Near Infrared Light (NIR) is a type of light just beyond the visible spectrum that can penetrate deeper into the skin than visible light.

The Nanoleaf LED Therapy Face Mask combines each light mode with NIR to stimulate cellular energy production and enhance blood circulation in the deeper layers of the skin for the most effective treatment and best results.

3

Are there any potential side effects to LED Therapy?

While light therapy is generally safe, in rare cases some individuals may experience side effects, including:

- Skin irritation
- Mild redness or swelling
- Increased sun sensitivity

If you have an existing skin condition or are on medications that increase light sensitivity, it's best to consult your doctor before starting LED Therapy. If you notice any unusual reactions during or after use, discontinue use immediately and seek professional advice.

4

How long do I use the mask for, and how often?

We recommend using the mask for 5-10 minutes at a time, 3-5 times per week.

- **5 minutes** offers a moderate dosage (30J/cm²), perfect for most skincare routines.
- **10 minutes** allows for deeper penetration at full intensity, enhancing results.

5

Can I use the mask more than the recommended time/amount?

Using the mask more than recommended will not speed up or improve your results, and in some cases might contribute to increased skin sensitivity. For the best results, stick to the suggested usage of 5-10 minutes, 3-5x per week.

6

How soon will I see results?

With consistent use, a noticeable improvement in skin texture and tone should be visible after 4 to 6 weeks. However, some users may begin to see initial improvements in as little as 2 weeks, with full benefits achieved after 8 to 10 weeks.

7

Can I use skincare products while using the mask?

Our mask is designed to be used on completely clean, dry, bare skin.

Applying skincare products before using the mask may create a barrier between the mask and the skin which could reduce the efficacy of the treatment. For this reason, we suggest applying your skincare products after each treatment session is complete.

8

Is this mask smart/app connected?

No, this mask is not a smart product or app-connected.

It features an attached controller that allows you to easily adjust treatment modes, intensity levels, and treatment time, giving you full control over your skincare routine without the need for a separate app.

9

What is the purpose of the eye attachments?

The eye attachments protect your eyes from direct light exposure during your treatment session to ensure a safe and comfortable experience.

10

Does the mask need to be plugged into power while I'm using it?

No! This mask is rechargeable and portable so you can move around freely without it needing to be connected to a power outlet during use. Just make sure to keep it charged by recharging it after every few sessions.

11

How do you clean the mask?

Simply wipe the mask with a soft damp cloth between sessions to keep it clean. If you wish, you can also use a small amount of mild detergent mixed with water on the cloth to remove any residual oils left behind on the mask by your skin. Do not put ANY component of the mask under running water or submerge it, or it will become damaged. (Make sure the device is always turned off and unplugged from power before cleaning)