6-in-1 Light Therapy Wand

User Manual

English

Français

Deutsch

Español


Italiano

6-in-1 Light Therapy Wand

User Manual

What's in the Box?


① 6-in-1 Light Therapy Wand

② Charging Cable

This device is rechargeable and portable. It can be used without being connected to an outlet.
Do not use the device while it is charging.

Controls

6 Treatment Modes

Anti-Aging

Red Light + Near Infrared (NIR) + Heating + Vibration
3 Intensity Levels

Acne

Blue Light + Vibration
3 Intensity Levels

Repair

Red Light + Cooling
3 Intensity Levels

Soothe

Heating + Vibration
3 Intensity Levels

Cold Therapy

Cooling at $15 \pm 3^{\circ}\text{C}$

Heat Therapy

Heating at $38 \pm 3^{\circ}\text{C}$

NOTE: It is normal for the wand's metal head to be warm after each treatment session.
Modes that utilize Cooling technology have an automatic 4 minute cooldown period between uses to avoid overheating.


Battery Status

Do not use the device while it is charging.

The 3 LEDs above the Intensity Button also help to indicate the battery level.


Low Battery

The 3 LEDs cycle from top to bottom, and beep every 5s.


Charging

The 3 LEDs cycle from bottom to top.


Fully Charged

All 3 LEDs stay on when plugged in and fully charged. Takes 4hrs for a full charge.


Important Things to Note


Warning: This is not a medical device and cannot be used for professional treatment or as a substitute for professional medical advice. The actual effect of use varies from person to person.


Avoid looking directly at the LED lights.

Do not drop or throw the product on hard surfaces.

Do not soak in water.


Clean with a damp cloth or with 70% alcohol.

How to Use

1. Wash your face before use.

Optional: Apply your favorite serums and moisturizers.

2. Hold the Power Button for 2s to Power On.

3. Press once more to select a treatment mode. The wand will beep at the 2 minute mark.

4. Once a treatment is done, the next mode will blink. Press the power button to start the next treatment mode, or select another mode.

5. Hold the Power Button for 2s to Power Off.

Suggested Wand Movements


Target specific areas of your face, or follow any of the suggested wand movements below.

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